# BALANSE LITE

## Chytrá váha

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**prove**

## SPECIFIKACE:

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Baterie: 2x AAA

Podpora: Bluetooth

LED displej: 75x31 mm

Rozměry: 280x280x22 mm

Barva: Modrá

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Hmotnost

Obezita

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Bílkoviny

Tělesný tuk

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Voda

Kostní hmota

Svaly

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**prove**

**SPECIFICATIONS:**

Material: ABS + tempered glass  
Maximum weight: 180 kg  
Measurement accuracy: ~0.1 kg  
Battery: 2xAAA  
Support: Bluetooth  
LED display: 75x31 mm  
Size: 280x280x22 mm  
Colour: Blue

**DEVICE OPERATION:**

1 measurement = 12 parameters

Weight

Obesity

Basal metabolic rate

Visceral fat

Body mass index

Protein

Body fat

Biological age

Water

Bone mass

Muscles

Live body weight

### 1. Installing the battery:

Remove the battery cover, insert 2 AAA batteries in the correct polarity.

### 2. Installing the app:

Scan the QR code above or download the FitGo app from the App Store (iOS) or Google Play (Android).

### 3. Registration:

Switch on Bluetooth on your mobile device, launch the FitGo app and follow the instructions in the app.

Up to 8 users can use one account on the device.

iOS App

Android App

## OPERATING INSTRUCTIONS:

1. Do not use the scale if you have implants in your body, such as a pacemaker, various kinds of prostheses, metal plates or screws. This can lead to malfunctioning of the device or inaccurate results. If in doubt, consult your doctor.
2. Stand flat on the surface of the scale without taking your feet off - otherwise, the device will only show the weight, without additional parameters.
3. Make sure the scale is on a flat, level and solid surface. Avoid weighing on carpet or other soft surfaces.
4. Do not let your feet or the surface of the scale get wet - it is very dangerous!
5. The scale is not protected against moisture, so keep it away from water.
6. Clean the scales with a slightly damp cloth - do not use chemical/abrasive cleaners.
7. If you do not plan to use the scale for a long time, remove the batteries.
8. The device is intended for home use only.

9. Children should use the scales under the close supervision of their parents or elders.
10. During the first weighing, the weight readings vary greatly.

The electronic scales will be automatically adjusted the first time you use them. The data of the first weighing includes the weight of the product itself, so the owner's exact data can be obtained when weighing again.

## RECYCLING RULES

The rules for recycling the stabiliser differ from those for household waste. Take the device to a specialised waste collection and recycling centre or contact the retailer from whom you purchased the device. Improper waste disposal is subject to penalties in accordance with the applicable laws of the country in which the product was purchased.

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